

# Australian Designed Safer Delivery System



## Thumb Grip Trocar

- The Patented "Thumb Grip Trocar" gives the Surgeon more control which helps prevent needlestick injuries during insertion.

## Superior Drainage

- Four flow channels for greater tissue contact. No holes means less tissue ingrowth and less risk of occlusion.

## Patient Comfort

- Silicone tubing allows for more flexibility
- The design means less tissue trauma on removal of drain
- Latex Free





## Features

- Radio opaque
- Available in round, flat and multiple channel
- With or without trocar
- All drains are double wrapped sterile, with a connector that will suit most types of wound drainage systems
- Compatible with high and low suction
- Visible skin mark






# Australian Designed, Latex Free Wound Drains

## Ordering Information Without Trocar

	Code	Description	Size	Case Qty
<b>Round Silicone</b>				
	VC500-32	Round Silicone Drain (with side holes)	10Fr	10
	VC500-33	Round Silicone Drain (with side holes)	16Fr	10
	VC500-34	Round Silicone Drain (with side holes)	20Fr	10
<b>Flat Silicone - Single Channel</b>				
	VC500-37	Flat Silicone Drain (with side holes) 3/4 length holes in drain section	7mm	10
	VC500-38	Flat Silicone Drain (with side holes) 3/4 length holes in drain section	10mm	10
<b>VariBlake™ Silicone - Multiple Channel</b>				
	VC500-41	VariBlake™ Silicone Drain without trocar (390mm channel section & 800mm tube without perforation)	10Fr	10
	VC500-42	VariBlake™ Silicone Drain without trocar (390mm channel section & 800mm tube without perforation)	15Fr	10
	VC500-43	VariBlake™ Silicone Drain without trocar (390mm channel section & 800mm tube without perforation)	19Fr	10
	VC500-46	VariBlake Silicone Drain without trocar (350mm channel section)	24Fr	10
<b>VariBlake™ Silicone - Multiple Channel Flat Drain</b>				
	VC500-44	Flat VariBlake™ Silicone Drain without trocar (390mm perforation & 800mm drain without perforation)	8mm	10
	VC500-45	Flat VariBlake™ Silicone Drain without trocar (390mm perforation & 800mm drain without perforation)	11mm	10

## Ordering Information With Trocar

	Code	Description	Size	Case Qty
<b>VariBlake™ Silicone - Multiple Channel</b>				
	VC500-51	VariBlake™ Silicone Drain with trocar (390mm channel section & 800mm tube without perforation)	10Fr	10
	VC500-52	VariBlake™ Silicone Drain with trocar (390mm channel section & 800mm tube without perforation)	15Fr	10
	VC500-53	VariBlake™ Silicone Drain with trocar (390mm channel section & 800mm tube without perforation)	19Fr	10
	VC500-53ST	VariBlake™ Silicone 19Fr Drain with 16Fr trocar (390mm channel section & 800mm tube without perforation)	19Fr	10
<b>VariBlake™ Silicone - Multiple Channel Flat Drain</b>				
	VC500-54	Flat VariBlake™ Silicone Drain with trocar (390mm perforation & 800mm drain without perforation)	8mm	10
	VC500-55	Flat VariBlake™ Silicone Drain with trocar (390mm perforation & 800mm drain without perforation)	11mm	10
<b>T Drain</b>				
	VC500-71	Silicone T Drain 10Fr drain only	10Fr	10
	VC500-72	Silicone T Drain 14Fr drain only	14Fr	10

™\* Trademark of International Medical Research.